
## **LUFRA**

# SAN MULTI VIOLET

A concentrated fragranced daily cleaner/ disinfectant suitable for use on all hard surfaces.

#### **DILUTION RATE**



#### **METHOD OF USE**



Fill an empty trigger bottle with 600ml of clean water.



Add 1 shot (30ml) of Lufra San Multi Violet to the trigger bottle.



Spray diluted Lufra San Multi Violet onto the area and work into the surface with a cloth.



For efficacy against bacteria, allow a 5 minute contact time on the surface.



Wipe away with a clean cloth. Leave the surface to air dry.

## EFFICACY STANDARDS

BACTERICIDAL - EN 1276			
CONTACT TIME	DILUTION	KILL RATE	
5 Minutes	1:23	99.999%	

ENVELOPED VIRUSES - EN 14476			
CONTACT TIME	DILUTION	KILL RATE	
60 Seconds	1:23	99.99%	

All testing methods conducted in dirty conditions by an independent laboratory.

Use biocides safely. Always read the label and product information before use.

#### **PRODUCT INFO**



Pack Size: 5L
Product Code: 810
Odour: Sweet
Colour: Yellow
pH (undiluted): 6-8
Shelf life: 2 years



For Health & Safety information, please scan the QR code provided.

THINK SAFE, CLEAN GREEN



### **LUFRA**

# SAN MULTI VIOLET

A concentrated fragranced daily cleaner/ disinfectant suitable for use on all hard surfaces.

#### **DILUTION RATE**



#### **METHOD OF USE**



Place a wet floor sign in the area to be cleaned.



Fill an empty bucket with 5 litres of clean water.



Dispense 6 pumps (180ml) of Lufra San Multi Violet into the bucket.



Immerse the clean mop head into the bucket of Lufra San Multi Violet solution and wring out the mop thoroughly.



Apply the mop to the floor and wipe in a figure of 8 motion, working the mop side to side.



Clean the dirty mop head in the bucket and repeat the previous steps until the floor is clean. Allow to air dry for a streak-free finish.

#### **EFFICACY STANDARDS**

BACTERICIDAL - EN 1276			
CONTACT TIME	DILUTION	KILL RATE	
5 Minutes	1:23	99.999%	

ENVELOPED VIRUSES - EN 14476			
CONTACT TIME	DILUTION	KILL RATE	
60 Seconds	1:23	99.99%	

All testing methods conducted in dirty conditions by an independent laboratory.

Use biocides safely. Always read the label and product information before use.

#### PRODUCT INFO



Pack Size: 5L
Product Code: 810
Odour: Sweet
Colour: Yellow
pH (undiluted): 6-8
Shelf life: 2 years



For Health & Safety information, please scan the QR code provided.

THINK SAFE, CLEAN GREEN

CHRISTEYNS
PROFESSIONAL HYGIENE



# LUFRA SAN MULTI VIOLET

#### **CLEANER/DISINFECTANT**



- KILLS BACTERIA
- EFFECTIVE AGAINST BACTERIA, YEAST AND ENVELOPED VIRUSES
- VIOLET FRAGRANCE

#### **PRODUCT DESCRIPTION**

LUFRA SAN MULTI VIOLET is a cleaner and disinfectant with a violet fragrance. LUFRA SAN MULTI VIOLET is neutral in pH with low foam and can be used by hand or through a scrubber drier.

Safe for use on all hard surfaces including stainless steel. Use to clean sinks, taps, baths, walls, doors, floors, fittings and fixtures, toilets, bidets and urinals. Suitable for disinfecting sports equipment, sunbeds and surfaces in changing rooms, saunas, steam rooms, etc.

#### **EFFICACY**

A concentrated cleaner-disinfectant. LUFRA SAN MULTI VIOLET has been tested against a range of standards, at a dilution as shown below by an independent laboratory.

STANDARD	DILUTION RATE	CONDITIONS	CONTACT TIME	PASSED
<u>Bactericidal</u>				
EN 1276	1:23 (4.2%)	Dirty	5 Minutes	≥5 log reduction
EN 13697	1:23 (4.2%)	Dirty	5 Minutes	≥5 log reduction
<u>Yeasticidal</u> EN 13697 EN 1650	1:23 (4.2%) 1:23 (4.2%)	Dirty Dirty	5 Minutes 5 Minutes	≥4 log reduction ≥3 log reduction
<u>Virucidal</u>				
<u>(enveloped virus)</u> EN 14476 (MVA)	1:23 (4.2%)	Dirty	60 Seconds	≥4 log reduction



#### **INSTRUCTIONS FOR USE**

Dilute with water. Mop, spray or scrub with LUFRA SAN MULTI VIOLET. Allow the contact time shown in the table and then rinse well with clean water.

Use biocides safely. Always read the label and product information before use.

#### **DILUTION RATE**

Cleaning, Deodorising & Sanitising (Bacteria and Yeast)	1:23 (4.2%)
Cleaning, Deodorising & Sanitising (Enveloped Viruses)	1:23 (4.2%)

#### **ADDITIONAL INFORMATION**

pH Level: 6-8

#### **STORAGE**

Use above 5°C. Store preferably above 0°C.

#### **PACKAGING**

Jerry Can 5 Litre

#### **NON WARRANTY**

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.

